



## Syringe holder



### Product description

Shielded Syringe Carriers reduces exposure while storing or transporting radioactive syringes. The overlapping lid design with snap-latch closure prevents streaming. There are two sizes and thicknesses of lead to choose from. The ends of the carriers are double thick to reduce the exposure from the ends of syringes.

### Specifications

#### **Type A: Shielded Syringe Carrier, Large**

*Dimensions: (LxWxH)*

- I.D: 21 x 7.6 x 7.4 cm
- O.D: 24 x 11.2 x 8.9 cm

*Lead Shielding:*

- Sides, top and bottom: 3mm
  - Ends: 6mm
- Weight: 5.1 kg

#### **Type B: Shielded Syringe Carrier, Small**

*Dimensions: (LxWxH)*

- I.D: 20.3 x 4.8 x 5 cm
- O.D: 23.5 x 8.6 x 6.6 cm

*Lead Shielding:*

- Sides, top and bottom: .3mm
  - Ends: .6,,
- Weight: 3.4 kg

#### **Type C: Shielded Syringe Carrier, Small**

*Dimensions: (LxWxH)*

- I.D: 19 x 4.3 x 4.3 cm
- O.D: 23.5 x 8.6 x 6.6 cm

*Lead Shielding:*

- Sides, top and bottom: 6mm
  - Ends: 13 mm
- Weight: 3.6 kg

#### **Type D: Shielded Syringe Carrier, Large**

*Dimensions:*

- I.D: 19.6 x 6.6 x 6.9 cm
- O.D: 24 x 11.2 x 8.9 cm

*Lead Shielding:*

- Sides, top and bottom: 6mm
  - Ends: 13 mm
- Weight: 7.7 kg

